powertube[®]

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Seite 3

Mode d'emploi

page 15

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pagina 27

Instruction

page 40

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
Instruction manual

This instruction manual is considered an integral part of the device. It shall be kept within reach at all times. Exact compliance with these instructions is essential for the correct handling of the device as well as for the safety of patient and user.


Tabel of Contents

		Page
1.	Illustration of the device	40
2.	Control panel with functions	41
3.	Description of the device	42
4.	Quick start guide	42
5.	Battery and battery compartment	43
6.	Specifications	44
7.	Handling	45
8.	Indications	46
9.	Contraindications	46
10.	Safety tips	47
11.	Care, maintenance and storage	47
12.	Errors and malfunctions, function test	47
13.	Accessories	48
14.	Guarantee	48
15.	Safety inspection	49
16.	CE-Conformity declaration	50

1. Illustration of the device

2. Control panel with functions

- 

Turning on the device (push button briefly)
- 

Select number of minutes (1-7 minutes)
- 

Start therapy (push button briefly) (manual STEP selection)
- 

Turning off the device (push button for an extended period of time or device will automatically turn off after 30 seconds)

3. Description of the device

Intended use

The POWERTUBE is intended for the electronic stimulation of nerves. It is used for self-treatment of symptoms and pain.

Its application does not replace treatment or instructions from a doctor, but has a positive influence on symptoms and pain.

Using the device can alleviate pain.

Function designation

The transcutaneous electronic nerve stimulation device POWERTUBE is an easy to use, battery-operated electronic device designed to stimulate the nerves.

Corresponding nerve points are stimulated directly with electronic impulses. The device can be operated by persons who do not have formal medical training.

The device is described in the enclosed easy to follow instruction manual to make self-treatment a straightforward process. The handling of the device is intended to suit users with no medical training. Various applications as well as corresponding therapy times are outlined in the instruction manual.

4. Quick start guide

Before using for the first time

Immediately after receipt, the device should be inspected for transport damage and be subject to a function test. In case of transport damage, the shipper shall be notified immediately and a damage report be drawn up to support any claim for damages. In addition to name and address of the recipient, the protocol shall include the date of receipt, type and serial number of the device as well as a description of the damage. The original packaging should be kept during the guarantee period to ensure that the device can be returned in its original packaging if necessary.


Every device is tested by the manufacturer in terms of functionality and safety before being shipped. Please observe the following to ensure that the device continues to operate properly following transport and installation:

The operator should only start using the device after carrying out a function test (see chapter 12).

5. Battery and battery compartment

Use alkaline batteries only and check the polarities (plus/minus) carefully. By using this kind of battery you may be certain to end your therapy time even in case the battery is almost empty by showing the sign «low battery». The POWERTUBE device is adjusted to this type of batteries.

If you use rechargeable accus (e.g. NiMH), the proper functioning of the sign «low battery» is not completely reliable because of the different mode of discharging. The POWERTUBE device may abruptly stop functioning without prior warning. If you still prefer to use rechargeable accus, please recharge them periodically.

6. Specifications

Display:	Time of therapy in minutes, 7 light-emitting diodes Therapy steps 1 - 3 Exchange light-emitting diode battery
Error messages:	Intermittent tone: handholds or contact areas are not applied properly Continuous tone: handholds have short circuited
Battery:	IEC 6LR61 9-volt block battery
Output voltage:	AC, maximum 9-volt pp and minimum 6-volt pp
Duration of treatment:	3 - 21 minutes
Programme sequence:	3 steps
Measurements:	Length: 18 cm, diameter: 3.1 cm
Weight:	≈ 1100 g
Application class:	BF 
Safety class:	Ila unit
Operation:	10°C - 40°C Humidity 30% - 75% without condensation, air pressure 700 hPa - 1060 hPa
Please refer to instruction manual!	
Production Date	month / year-serial number

6. Handling

1. Turning on the device: push I/GO button briefly.
2. Set appropriate therapy time per step by means of TIME button. Start with one minute per step, i.e. the complete treatment time is 3 minutes. Increase the time with every treatment by 1 minute per step. Complete steps 1 and 2 (see below and beside) once a day.
3. Place device on trigger point.
4. Start treatment by pushing the I/GO button.
5. You will hear a beep signal after completion of every step.
6. The end of the treatment is announced with two beep signals.
7. Turning off the device: push I/GO button for a prolonged period of time or device turns off automatically after 30 seconds.

Error message: Insufficient or no skin contact is indicated with continuous beep signals.

Step 1 Local treatment at trigger point

Step 2

Treatment of abdominal area

Please ensure good body contact. Moisten skin with water if necessary.

Consume large quantities of water before and after the treatment.

8. Indications

- Pain relief

9. Contraindications

- Pregnancy: do not use the device
- Patients with a cardiac pacemaker: do not use the device
- Open wounds: do not use the device
- Fresh scars: do not use the device
- Do not use in combination with other devices

10. Safety tips

- POWERTUBE may only be used for the electronic stimulation of nerves and musculature
- Do not use the device for patients who have been fitted with a cardiac pacemaker
- Simultaneous use with other devices is prohibited
- POWERTUBE may only be operated with a 9-volt battery
- Remove the battery if the device is not being used for a prolonged period of time
- POWERTUBE may only be operated with the accessories included in the range of delivery
- Famedtec GmbH is only responsible for the reliability and safety of POWERTUBE if:
- POWERTUBE is operated in accordance with the instruction manual
- expansions and modifications are carried out only by persons who are authorised by Famedtec GmbH
- Please recycle the batteries in accordance with national regulations
- POWERTUBE should be disposed of as hazardous waste

11. Care, maintenance and storage

- Effective protection of the POWERTUBE from damage includes:
- Appropriate operation and maintenance
- Protection from moisture and dirt
- Wipe electrodes after every use to disinfect them
- POWERTUBE should be stored in a dry place

12. Errors and malfunctions, function test

- Discharge of the battery is displayed on a battery indicator LED

- After being turned on, POWERTUBE automatically runs a battery test
- An acoustic warning signal during therapy indicates that the electrodes have not been properly applied

Test the proper functioning of your POWERTUBE device yourself

1. Insert a new or newly charged battery into the battery compartment.
2. Connect the handheld with the cable for connection to the POWERTUBE device.
3. Push the start-button on the control panel.
4. Let the standard setting of 1 minute therapy time stand.
5. Start the therapy by pushing the start-button again and touch with the handheld the POWERTUBE casing (not the therapy head!).
6. If the POWERTUBE device now emits a continuing peep-sound, the device functions perfectly.
7. Finish the test by pushing the start-button until the device is turned off completely.

13. Accessories

- Hand electrodes
- Instruction manual
- IEC 6LR61 9-volt battery
- Portable case for POWERTUBE
- Handle ball gold
- Handle spoon gold

14. Guarantee

- Device and accessories shall be inspected for transport damage immediately after receipt. Corresponding claims for damages can only be asserted if the vendor or shipper have been notified immediately. A damage report shall be compiled.
- The guarantee for POWERTUBE is two years from the day of delivery.

- The extent of guarantee includes free repair of the device, provided the damage was caused by a material or manufacturing error.
- Additional claims, particularly claims for damages, are excluded.
- The repair may only be carried out by Famedtec GmbH, one of its agencies or authorised specialist dealers. The guarantee shall become void if improper modifications or repairs are carried out.
- The guarantee period will neither be extended nor renewed because of guarantee work.
- Famedtec GmbH reserves all rights to this instruction manual, particularly the right of duplication and distribution as well as translation. No part of this instruction manual may be copied or distributed in any form without prior written consent from Famedtec GmbH.
- Information contained in this instruction manual may be modified or upgraded by Famedtec GmbH.

15. Safety inspection

- § 6 of the directive on installing, operating and using active medical products stipulates regular safety inspections to guarantee the safety of the device against age, wear and tear etc.
- The operator is responsible for having these safety inspections carried out on the device to the prescribed extent and in due time. Safety inspections may only be carried out by the manufacturer or by explicitly authorised persons.
- The following safety inspections are required for POWERTUBE:
- Examination of labelling and instruction manual
- Visual inspection of device and accessories for damages
- Function test of all sensors and control lamps on the device

16. EG-Conformity declaration

CE 0482

Name and address of the company	Famedtec GmbH Breiten 9, CH-3636 Forst
We declare in the position as solely responsible company that	
the medical electrical equipments	POWERTUBE (gold) Typ FPT.B1 POWERTUBE (silver) Typ FPT-D.B2
in safety class	Ila according to Appendix IX of the Directives 93/43/EWG
all required standards of the directives mentioned in the medical dispositions 93/42/CEE are herewith completely fulfilled and adhered to.	
Conformity evaluation procedure	Medical dispositions 93/42/CEE, Appendix II
Responsible evaluation centre for the issuance of this conformity declaration	MEDCERT GmbH Pilatuspool 2 D-20355 Hamburg CE 0482
This conformity declaration is valid only together with the corresponding clearance documents for the production of theses devices.	
CH-Forst BE 2020	Werner Frischknecht (CEO of Famedtec GmbH)
Valid time of the certificate	2023-01-23

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Instruction Manual
Health Canada

Instruction Manual
TGA Australia

© 2020 Famedtec GmbH. All rights reserved. All trade names are registered trademarks of

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